B4®

Pre-Impression Surface Optimizer

DIRECTIONS FOR USE – ENGLISH

For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

B4® Pre-Impression Surface Optimizer is designed to be used as a pre-impression conditioning solution. When applied to substrates such as sulcular tissue, tooth preparations, ceramic, composite or metal core build-up material or implant abutments or copings, it provides a uniform surface condition with reduced surface tension prior to the application of all Aquasil® brand Smart Wetting® or Tissue Managing Impression Materials. This standardized surface with reduced contact angle on all impressed substrates allows capture of fine preparation detail both supra and subgingivally.

1.1 INDICATIONS

Preparation conditioning solution for uniform surface wetting prior to elastomeric impression procedures.

1.2 CONTRAINDICATIONS

B4® Optimizer is contraindicated for direct application to dental pulp tissue.

1.3 DELIVERY FORMS

(Some delivery forms may not be available in all countries)

B4® Optimizer is available in:

- 6mL tri-laminated bottle

1.4 COMPOSITION

Polyethylene Glycol (PEG); Ethanol; Water; Surfactant; Pigments

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 WARNINGS

1. B4® Optimizer contains PEG (polyethylene glycol). Avoid prolonged or repeated contact with skin (allergic contact dermatitis), oral soft tissues, and eyes. Do not take internally.

Eye contact: B4® Optimizer contains PEG which may be irritating to eyes. Before using this product wear protective glasses as well as covering the patient's eyes to protect from splashing material. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.

Skin contact: B4® Optimizer contains PEG which can cause skin sensitization (allergic contact dermatitis) in susceptible individuals. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water after contact. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention.

Oral mucosa contact: Avoid prolonged contact with oral soft tissues. If accidental contact occurs, flush mucosa with plenty of water and expectorate water. If sensitization of mucosa persists, seek medical attention immediately.

2. B4® Optimizer should not be used with patients who have a history of severe allergic reaction to any of the components.

2.2 PRECAUTIONS

- This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- The bottles cannot be reprocessed. To prevent bottles from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves. During treatment clinicians with patient contact should not handle bottles.
- The B4® Optimizer bottle should be tightly closed immediately after use.
- Use only in well ventilated areas.
- Flammable: B4® Optimizer contains ethanol. Keep away from sources of ignition.

Interactions:

- B4® Optimizer should not be used in conjunction with adhesive procedures. Marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. Rinse areas thoroughly following use.

2.3 ADVERSE REACTIONS

1. Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen

on the skin. Mucous Membranes: inflammation, edema, sloughing. (See Warnings)

2. B4® Optimizer may cause pulpal effects. (See Contraindications)

2.4 STORAGE

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Keep out of direct sunlight and store in a well ventilated place at room temperature between 10°-24°C/50°-75°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS FOR USE

PRE-IMPRESSION PREPARATION CONDITIONER

- 3.1.Preparation: Complete tooth preparation, including core buildup or implant abutment/coping connection and/or retraction cord placement as necessary, in usual fashion.
- 3.2.Cleaning: Remove provisional restoration and residual temporary cement (if applicable). Clean preparation with a rubber cup and pumice or a cleaning paste such as Nupro® Prophylaxis Paste. Wash thoroughly with water spray and air dry. Clean freshly instrumented preparation or implant component with water spray and then air dry.
- 3.3.Cord placement: If retraction cord is to be used, place per usual technique in accordance with manufacturer's instructions.
- 3.4.Application of B4® Optimizer:
 - 3.4.1Dispense 2-3 drops of B4® Optimizer from the bottle into a clean dappen dish. Apply additional drops for multiple preparations, if needed. Replace cap promptly.
 - 3.4.2 Using a clean disposable brush, apply B4® Optimizer to the entire preparation/abutment/coping and sulcus, thoroughly (re)wetting the packed retraction cord, if applicable. Air-thin for a uniform film but do not rinse or dry applied B4® Optimizer.
 - 3.4.3Remove cord, if applicable, and confirm hemostasis per usual technique. If additional chemical agents are used to establish hemostasis, thoroughly rinse preparation and sulcus and re-apply B4® Optimizer as outlined in section 3.4.2 once bleeding is controlled. Once hemostasis is confirmed, proceed immediately to impression as per Aquasil® brand Impression Material's Directions for Use.
- 3.5.After impression completion, thoroughly rinse areas where B4 Optimizer was applied.

4. HYGIENE

Cross-contamination - Infection

Do not reuse single use products. Dispose of in accordance with local regulations.

4.1 BOTTLES

Cross-contamination - Infection

- The bottles cannot be reprocessed. Dispose of contaminated bottles in accordance with local regulations
- To prevent bottles from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves. Do not reuse bottles if contaminated.
- During treatment clinicians with patient contact should not handle the bottle

Incidental contact of the bottle with water, soap or a water-based hospital-level disinfection solution will not damage bottle body. Do not allow any solution contact with contained material. Discard material that has been in contact with any fluid or non-sterile instrument. Repeated liquid contact may damage label. Dry the bottle with a lint-free single-use cloth.

NOTE: Vigorous wiping can destroy the label. Wipe bottles gently.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard is used: "YYYY-MM-DD"
2. The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot number on bottle
 - Expiration date on bottle

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