NRC™

Non-Rinse Conditioner

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Conditionneur sans rinçage

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Grabado ácido sin lavado

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CAUTION: For dental use only.

USA: Rx only.

Content	Page
1 Product description	2
2 Safety notes	3
3 Step-by-step instructions	4
4 Hygiene	5
5 Lot number and expiration date	5

1 Product description

NRC™ is a self-priming conditioning liquid used in combination with Prime&Bond® NT Nano-Technology Dental Adhesive for Dyract® compomer restorations where enhanced adhesion to enamel is required.

NRC™ offers a fast way to condition and prime enamel and dentin in an easy application procedure without the need to rinse with water.

1.1 Delivery forms

- 4.5 ml bottle.

1.2 Composition

- Maleic acid
- Itaconic acid
- Water

1.3 Indications

NRC™ is to be used prior to application of Prime&Bond® NT and Dyract® compomer materials in the following situations:

- Stress-bearing Class I and II permanent Dyract® compomer restorations.
- Class IV Dyract® compomer restorations.
- Cavities where beveled enamel margins are preferred for esthetic reasons.

NRC™ may be used for all other Dyract® compomer applications, if wished.

1.4 Contraindications

- Application to dentin close to the pulp (less than 1 mm).
- Use in combination with adhesives other than Prime&Bond® NT.

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

The material contains organic acids which may cause burns and may be irritating to skin, eyes and oral mucosa. The material may cause sensitization in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation. If contact with skin occurs, remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use. Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Tightly close bottle immediately after use.
- Use in a well-ventilated area. Avoid inhaling spray.
- To prevent the bottle from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves.

2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Sensitization. Irritation.
- Contact with mucous membranes: Inflammation (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Keep the product out of direct sunlight and store in a well-ventilated place.
- Store at temperatures between 10 °C and 24 °C (50 °F and 75 °F). Use the product at room temperature.
- Protect from moisture.
- Do not freeze.
- Do not use after expiration date.

3 Step-by-step instructions

3.1 Cleaning

Surface cleanliness is paramount for the development of adhesion.

1. Clean uninstrumented enamel and dentin with a rubber cup and pumice or a non-fluoride cleaning paste such as Nupro® prophylaxis paste.
2. Preparing a fresh surface with a finishing bur will significantly increase bond strength to enamel.
3. Wash thoroughly with water spray and then air-dry.
4. Remove rinsing water completely by blowing gently with an air syringe or by blot drying with a cotton pellet.
5. Do not desiccate the tooth structure.

3.2 Pulp protection

1. In deep cavities cover the dentin close to the pulp (less than 1 mm) with a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see complete Instructions for Use) leaving the rest of the cavity surface free for bonding.

3.3 Application

1. Dispense NRC™ into a Dentsply Sirona Applicator Dish¹ or standard dappen dish.
2. Apply sufficient amounts of the material to enamel and dentin with an Applicator Tip¹ or disposable brush.
3. Leave undisturbed for 20 seconds. Do not rinse.
4. Remove excess material by blowing gently with an air syringe or blot drying with a cotton pellet. Do not desiccate the tooth structure.

Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat the application of NRC™.

3.4 Application of Prime&Bond® NT

1. Dispense Prime&Bond® NT adhesive directly onto a fresh Applicator Tip or onto a disposable brush. Alternatively, dispense into a fresh Dentsply Sirona Applicator Dish or standard dappen dish.
2. Immediately apply ample amounts of Prime&Bond® NT adhesive to thoroughly wet all cavity surfaces. These surfaces should remain fully wet for 20 seconds and may necessitate additional application of mixed adhesive.
3. Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance. If not, repeat application and air-dry.
4. Light cure for a minimum of 10 seconds². Ensure uniform exposure of all cavity surfaces.
5. Immediately place restoration with luting compomer over the cured Prime&Bond® NT adhesive.

¹ Dentsply Sirona Applicator Dish and Applicator Tips are available from your dental dealer.

² Check curing light for minimum curing output of at least 300 mW/cm².

4 Hygiene

Cross-contamination.

Infection.

1. Do not reuse single use products. Dispose in accordance with local regulations.

5 Lot number () and expiration date ()

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date

© Dentsply Sirona 2014-12-05

[These Instructions for Use are based on Master Version 01]

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Dentsply DeTrey GmbH