core-X[®] flow

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core-X[®] flow

Dual Cure Core Build-Up Material and Cement for Endodontic Posts

CAUTION: For dental use only.
In the USA: Rx only

Content	Page
1 Product description	2
2 General safety notes	3
3 Step-by-step instructions	4
4 Hygiene	7
5 Lot number and expiration date	7

1 Product description

core-X[®] flow consists of two-components, base and catalyst, which when mixed form a dual cured, highly filled, composite resin core build-up and post-cementation material.

core-X[®] flow uses a biocompatible urethane resin and is supplied in a tooth-colored shade which is ideal in situations where esthetics and show-through of the core are of primary concern.

1.1 Delivery forms

core-X[®] flow is available in a convenient dual-barreled syringe. This delivery system simplifies direct, precise intraoral delivery and minimizes product waste.

1.2 Composition

- Urethane Dimethacrylate
- Di- & Tri-functional Methacrylates
- Barium Boron Fluoroaluminosilicate Glass
- Camphorquinone (CQ) Photoinitiator
- Photoaccelerators
- Silicon Dioxide
- Benzoyl Peroxide

1.3 Indications

- Vital or non-vital tooth core build-up (replacement of existing restorations and/or lost tooth structure) as a base prior to fabricating an indirect restoration.
- Cementation of endodontic fiber-posts.

1.4 Contraindication

- **core-X[®] flow** is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

1.5 Compatible adhesives

core-X® flow is to be used after the application of a suitable dentin and enamel adhesive.

core-X® flow is compatible with numerous dentin/enamel adhesive systems including all Dentsply Sirona adhesives designed for use with dual/self cured composite restoratives (see complete Directions for Use of selected adhesive). Other dentin/enamel adhesive systems may be used at the discretion and sole responsibility of the dental practitioner.

2 General safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Directions for Use.

2.1 Warnings

- **core-X® flow** contains polymerizable methacrylates which may irritate to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitisation or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitisation of mucosa persists, seek medical attention.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
3. Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
4. Syringe should be tightly closed by replacing the original cap immediately after use.
5. **core-X® flow** is intended to be used as a dual cure product. To obtain optimal curing follow the step-by-step instructions. Self cured, dual cured and light cure methods are integrated into the procedural steps.
6. **core-X® flow** should be mixed in equal volumes directly through the supplied mix tips for optimum performance. Variations from this 1:1 ratio may adversely affect working time and compromise final physical properties. Dispensing and mixing other than through the supplied mixing tips is not recommended.
7. **core-X® flow** should extrude easily. **DO NOT USE EXCESSIVE FORCE.** Excessive pressure may result in unanticipated extrusion of the material or cause syringe rupture.
8. **core-X® flow** when used for post-cementation has a working time of only 40 seconds! Prepare and organize the necessary materials and instruments to allow a fast insertion of the post.

9. Interactions

- Eugenol-containing dental materials should not be used with **core-X® flow** because they may interfere with hardening and cause softening of the polymeric components of the materials.
- Variable in-vitro data exist regarding use of light cured-only adhesives without Self Cure Activator in conjunction with dual cured or self cured composite restorative materials such as **core-X® flow** in limited or no light curing applications. Chemical/product incompatibility may adversely affect product efficacy, leading to premature restoration failure.
- As with any dual cure resin system, the use of a dual cure adhesive system can shorten working time. This effect should be investigated in the laboratory prior to clinical use.
- The **core-X® flow** is a light cured material. Therefore, it should be protected from ambient light. Proceed immediately once material has been placed.

2.3 Storage conditions

core-X® flow syringes should be stored with the original cap only, at temperatures between 10 °C/50 °F and 24 °C/75 °F. Keep out of direct sunlight and protect from moisture. Refrigerated storage is acceptable when not in use. Allow material to reach room temperature prior to use. Do not freeze. Do not use after expiration date.

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

2.4 Adverse reactions

Product may irritate the eyes and skin. **Eye contact:** Irritation and possible corneal damage. **Skin contact:** Irritation or possible allergic response. Reddish rashes may be seen on the skin. **Mucous Membranes:** Inflammation (see 2.1 Warnings).

3 Step-by-step instructions

Working and Setting Times for Core Build-Up and Endodontic Post Cementation

Working Time at body temperature (37 °C) – for manipulation/post seating after intraoral placement	40 seconds
Working Time at room temperature (22 °C) – for application after start of mix	1 minute 30 seconds
Setting Time Light Cure	20 seconds ¹ (3 mm depth)
Setting Time Self Cure (catalyst initiated)	2-3 minutes

3.1. Core Build-Up

3.1.1 Cleaning and Preparation of Tooth

1. Finish preparation by removing existing restorations and/or caries.
2. Clean freshly cut enamel and dentin with water spray and air dry. Do not desiccate dentin.
3. Place pin, post or matrix as needed.

3.1.2 Pulp protection, tooth conditioning/dentin pretreatment, adhesive application

Refer to applicable manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of **core-X® flow**.

3.1.3 Dispensing and mixing core-X® flow dual-barreled syringe

Danger of injury due to excessive force.

1. Apply slow and steady pressure on the syringe plunger.
2. Do not use excessive force – syringe rupture may result.

1. If refrigerated prior to use, remove from refrigerator and bring to room temperature.
2. Remove syringe cap and dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross-contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
3. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mix tip with the v-shape notch on the syringe flange. Turn colored mix tip cap 90 degrees in a clockwise direction to lock in place on syringe.
4. Gently depress syringe plungers to begin the flow of material. **DO NOT USE EXCESSIVE FORCE.** If force is encountered, remove syringe from operating field, remove and discard mix tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mix tip as outlined above. Dispense a small amount through the mix tip onto a mixing pad and discard.
5. Proceed immediately to placement, **Step 3.1.4. Placement of core-X® flow.**

3.1.4 Placement of core-X® flow

core-X® flow produces a syringable mixed consistency for direct placement using the dual-barreled syringe and a banded matrix retainer or clear crown form. Alternatively, mixed material may be expressed onto a clean paper pad and manually placed with the operator's placement instrument(s) of choice. Working time of mixed material at room temperature (22 °C/72 °F) is approximately 1'30" (one minute thirty seconds) when protected from ambient light. Elevated oral temperatures and/or contact with dual cure adhesives may shorten work time (**see 2.2 Precautions**).

Any of the following placement techniques are acceptable. The preferred method is to use both visible light and chemical self cure (**see Step 3.1.4.3 Dual Cure**).

1. Visible Light Cure

Place mixed material directly in 2-3 mm increments and light cure increment for 20 seconds¹.

2. Chemical Self Cure

Place large increment(s) or load mixed material into a clear crown form and seat onto the preparation. Allow to self cure approximately 2-3 minutes. Visible light curing¹ of external surface is optional, but advisable.

3. Dual Cure

After placing the first 2-3 mm increment, and visible light curing for 20 seconds¹, large increments may then be placed or material loaded into a clear crown form, seated onto the preparation and allowed to self cure. Visible light curing¹ of external surface is optional, but advisable.

3.1.5 Finishing the core-X® flow

After initial self cure set and removal of matrix, a surface light cure of the material for an additional 20 seconds¹ is recommended. Initial crown preparation (gross reduction) may begin immediately using rotary instruments. If used per **step 3.1.4.2 (Chemical Self Cure)** or **3.1.4.3 (Dual Cure)**, utilizing chemical self cure, greatest bulk hardness is obtained 10 minutes after start of mixing at which time the final preparation may be completed.

3.2 Endodontic Post Cementation

Tooth Conditioner/Dentin Pretreatment/Adhesive application.

Refer to applicable manufacturer's directions for tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of **core-X® flow**, as outlined below.

Note: If core build-up will immediately follow post cementation, treat all available tooth structure beyond post space with adhesive for core build-up bonding.

1. Preparation of post:

Treat surface of post according to adhesive manufacturer's and/or dental laboratory's instructions, i.e., adhesive application.

2. Dispensing and mixing **core-X® flow** dual-barreled syringe.

1. Prepare syringe, bleed material and install mix tip as outlined above **3.1.3. (Dispensing and mixing core-X® flow dual-barreled syringe 1-4).**

3. Proceed immediately with placement. Dispense the same material directly from the syringe onto the post surface and to the orifice of the post hole preparation.

4. Immediately place post into final position and stabilize.

Important note: Work time of applied cement within a post space pretreated with dual cure adhesive is 40 seconds maximum from application. Post must be fully seated and stabilized within 40" (**see 2.2 Precautions, 9. Interactions**).

Inadequate insertion of post due to short working time.

1. At body temperature remaining working time is 40 seconds from start of application.
2. Prepare and organize necessary material and instruments to allow a fast insertion of the post.

Once stabilized, light cure, using a visible light curing unit designed to cure CQ initiated methacrylates (spectral output including 470 nm), with a minimum output of 550 mW/cm². Cure all accessible areas of the post for 20 seconds.

Inadequate polymerization due to insufficient curing.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output before each procedure.

5. Proceed with the core build-up as described in **chapter 3.1 (Core Build-Up)**.

4 Hygiene

4.1 Cleaning the syringe

For dual-barreled syringe, remove used mix tip and discard appropriately. Replace original syringe cap prior to storing. The dual-barreled syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.

4.2 Disinfecting the syringe

Disinfect the syringe with a hospital-level, tuberculocidal disinfectant solution according to national/local regulations. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Some phenolic-based agents and iodophor-based products may cause surface staining. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

5 Lot number () and expiration date ()

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date

© Dentsply Sirona 2007-11-28

[These Directions for Use are based on Master Version 01]

¹ Use suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 550 mW/cm² exposure. Refer to curing light manufacturer's recommendations for compatibility and curing recommendations.

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Dentsply DeTrey GmbH